

血细胞分析仪用校准品 Hematology Calibrator

校准品靶值表 Zybio Calibrator target value and range

LOT H201201
2021-02-02

仪器型号 Instrument Model	参数 Parameter	靶值 Target Value	不确定度 Acceptable Range
Z3、Z3 CRP	WBC $\times 10^9/L$	8.90	± 0.26
	RBC $\times 10^{12}/L$	4.55	± 0.11
Z30、Z30 CRP	HGB g/L	144	± 2
Z31、Z31 CRP	MCV fL	89.8	± 2.2
	PLT $\times 10^9/L$	265	± 21
Z5、Z50、Z51、Z52	WBC $\times 10^9/L$	9.05	± 0.28
	RBC $\times 10^{12}/L$	4.57	± 0.11
	HGB g/L	138	± 2
	MCV fL	85.5	± 2.2
	PLT $\times 10^9/L$	260	± 21

注意事项：在使用血细胞分析仪用校准品之前，请阅读包装说明书。如果未严格遵循说明，可能会导致校准错误。
血细胞分析仪用校准品开瓶有效期为7天。

Notes : Before using Hematology Calibrator, read the insert. If it is not strictly followed, calibration errors may be caused.
Once opened, the Zybio Hematology Calibrator are stable for 7 days.



Zybio Inc.

住所：重庆市大渡口区建桥工业园C区太康路6号30栋第1-4层 邮编：400082 联系方式：023-68190901
E-mail: zy@zybio.com 网址: http://www.zybio.com
Add: Floor 1 to Floor 4, Building 30, No.6 of Taikang Road, Block C of Jianqiao Industrial Park, Dadukou District, Chongqing, China 400082 Tel:+86(0)23 6865 5509 E-mail: info@zybio.com Web: www.zybio.com