

CBC-CAL PLUS

For Zybio
HEMATOLOGY CALIBRATOR

CALIBRATOR

LOT PLUS0121



2021-02-05

Assay values and ranges provided by Zybio

Instrument Model	Parameter	Target Value	Acceptable Range
Z3、Z3 CRP	WBC x10 ⁹ /L	8.90	± 0.26
	RBC x10 ¹² /L	4.55	± 0.11
Z30、Z30 CRP	HGB g/L	144	± 2
Z31、Z31 CRP	MCV fL	89.8	± 2.2
	PLT x10 ⁹ /L	265	± 21
Z5、Z50、Z51、Z52	WBC x10 ⁹ /L	9.05	± 0.28
	RBC x10 ¹² /L	4.57	± 0.11
	HGB g/L	138	± 2
	MCV fL	85.5	± 2.2
	PLT x10 ⁹ /L	260	± 21

Notes: Before using CBC-CAL PLUS, read the package insert. Calibration errors may result if instructions are not followed exactly.
CBC-CAL PLUS for Zybio has 7-day open vial stability



Zybio Inc.

Add: Floor 1 to Floor 4, Building 30, No.6 of Taikang Road, Block C of Jianqiao Industrial Park, Dadukou District, Chongqing, China 400082 Tel: +86-23-6869 0698 Email: info@zybio.com Web: www.zybio.com