

# 血细胞分析仪用校准品 Hematology Calibrator

## 校准品靶值表 Zybio Calibrator target value and range

LOT 510101  
 2021-02-10

| 仪器型号<br>Instrument Model | 参数<br>Parameter        | 靶值<br>Target Value | 不确定度<br>Acceptable Range |
|--------------------------|------------------------|--------------------|--------------------------|
| Z3、Z3 CRP                | WBC $\times 10^9/L$    | 8.90               | $\pm 0.26$               |
|                          | RBC $\times 10^{12}/L$ | 4.55               | $\pm 0.11$               |
| Z30、Z30 CRP              | HGB g/L                | 144                | $\pm 2$                  |
| Z31、Z31 CRP              | MCV fL                 | 89.8               | $\pm 2.2$                |
|                          | PLT $\times 10^9/L$    | 265                | $\pm 21$                 |
| Z5、Z50、Z51、Z52           | WBC $\times 10^9/L$    | 9.05               | $\pm 0.28$               |
|                          | RBC $\times 10^{12}/L$ | 4.57               | $\pm 0.11$               |
|                          | HGB g/L                | 138                | $\pm 2$                  |
|                          | MCV fL                 | 85.5               | $\pm 2.2$                |
|                          | PLT $\times 10^9/L$    | 260                | $\pm 21$                 |

**注意事项：**在使用血细胞分析仪用校准品之前，请阅读包装说明书。如果未严格遵循说明，可能会导致校准错误。  
血细胞分析仪用校准品开瓶有效期为7天。

**Notes :** Before using Hematology Calibrator, read the insert. If it is not strictly followed, calibration errors may be caused.  
Once opened, the Zybio Hematology Calibrator are stable for 7 days.

## Zybio Inc.

住所：重庆市大渡口区建桥工业园C区太康路6号30栋第1-4层 邮编：400082 联系方式：023-68190901  
E-mail: zy@zybio.com 网址: <http://www.zybio.com>  
Add: Floor 1 to Floor 4, Building 30, No.6 of Taikang Road, Block C of Jianqiao Industrial Park, Dadukou District,  
Chongqing, China 400082 Tel:+86(0)23 6865 5509 E-mail: info@zybio.com Web: www.zybio.com

